

**The choice
is yours!**

 **atair**®



The choice is yours!

StarMedTec presents **atair®**, its first high power diode laser at a wavelength of 980 nm. A maximum power output of 100 W perfectly suits this laser system for a fast, efficient coagulation of large tissue volumes (e.g. kondyloma, tumors). Thermal behaviour and haemostasis are significantly improved compared to Nd:YAG.

The **atair®** is based on StarMedTec's new platform of laser systems with all its innovative features: Ergonomic and compact design, reliable and state-of-the-art components, automatic fiber recognition and an intelligent, intuitive user interface on a brilliant color touchscreen.

Air cooled all-solid-state-technology ensures smooth operation and minimal maintenance.

A broad range of applicators is available by StarMedTec.



Powerful,
reliable
technology

Intelligent,
intuitive
user interface

Compact,
ergonomic design



**Quick solutions, clear answers:
the StarMedTec hotline is there for you.**

+ 49 (0) 8151/268 61-0

StarMedTec GmbH

Kreuzstraße 22 | 82319 Starnberg

Germany

info@starmedtec.com

www.starmedtec.com